

<b>FLUOCYNE® 10%</b> <b>Sodium fluorescein</b>
---

**500 mg per 5 mL injectable ampoule**

**Box of 10 ampoules**

FLUOCYNE® 10% is used by intravenous injection for the examination of retinal superficial vasculature (the examination of retinal deep vasculature is performed with Indocyanine Green).

FLUOCYNE® 10% (Sodium fluorescein) is used in the diagnosis of retinal pathologies, vascular pathologies, diabetic retinopathy.

Sodium fluorescein is also used in the examination of the terminal vasculature of the tissues. In the diagnosis of vascular pathologies, Sodium fluorescein facilitates the differentiation of healthy from diseased or damaged tissue.

10% Sodium fluorescein can be diluted in sterile water for injectable preparations and be used as a 2% solution for detecting corneal lesions and foreign bodies. 5% Sodium fluorescein is given intravenously for treating acrocyanosis and Raynaud's disease.

**Dosage: 200 to 500 mg by intravenous injection diluted or not.**